


No. 08/PV/GMP
Section : GMP
DT. 16.03.2021

Sub: Providing data for SDOTE/OTE cases in vendor registration portal.

Ref: GMs approval No. 3210/Gen/PV/M1 dated 12/02/2021

With reference to the above, the requisite information related to manufacturing of technology and testing/inspection facility of Mat-O-Bond are give below:

Sl. No.	Manufacturing Technologies Required	
01	Vacuum Distillation	
	Name of Item of Manufacturing	
01	Round bottom Flask 4 neck with cap set (20 Ltr.)	
02	Methanol Received Conical Flask	
03	Thermometer	
04	Teflon Rod for Stirrer	
05	Teflon cap for Bottom Flask	
06	Hot Water Heater set	
07	Plastic cane Air tight	
08	Sample container	
09	Nitrogen Cylinder gauge/Valve	
10	Nitrogen Cylinder UHP	
11	Weighing container for MAPO	
12	Condenser (Water Flowing Tube)	
13	Deep Freezer for storage of MOB	
	Testing Apparatus/Facility	
01	Volatile Mater (Oven Hot air 100°C with temp. Balance)	
02	Moisture Content (Karl Fisher Aqua meter)	
03	OH Value (Itrater) Heating Mental Plate	


Sr. DO/GMP

DO/PV

Copy:

DO/WI : for information please.

DO/Lab: for information please.